

# Programmier-Verbindungen auf einen Blick

	NoahLink Wireless	ConnexxAir	Flexstrip-Adapter	Prog.adapter 10	Prog.adapter 312	Prog.adapter 13
Insio Nx IIC	-	✓	-	-	-	-
Insio Nx CIC	-	✓	✓	-	-	-
Silk Nx	-	✓	✓	-	-	-
CROS Silk Nx	-	✓	✓	-	-	-
Pure 10 Nx	-	-	-	✓	-	-
Pure 312 Nx	✓	-	-	-	✓	-
Pure 13 Nx	✓	-	-	-	-	✓
Pure Charge&Go Nx	✓	-	-	-	-	-
Motion 13 Nx	✓	-	-	-	-	✓
Motion 13 P Nx	✓	-	-	-	-	✓