

# Programmierverbindungen auf einen Blick

Überblick über die Kompatibilität unserer Programmierverbindungen mit unseren Hörsystemen.



		Styletto Connect	Styletto	Pure 10 Nx	Pure 312 Nx	Pure 13 Nx	Pure Charge&Go Nx	Motion 131 Nx	Motion 13P Nx	Motion Charge&Go Nx	Silk Nx	Insio Nx (IC)	Insio Nx (CIC)	Insio Nx (TC/TE)	Ace primax	Pure primax	Cellion primax	Carat primax	Motion SP primax	Motion P primax	Motion SX primax	Motion SA primax	Motion S primax	Silk primax	Insio primax	Intuis 3 RIC 312	Intuis 3 S	Intuis 3 M	Intuis 3 P/SP	Intuis 3 IdO	Intuis 3 Click ITE	
KABELGEBUNDENE PROGRAMMIERADAPTER	Programmier-Adapter 10	-	-	•	-	-	-	-	-	-	-	-	-	-	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
	Programmier-Adapter 312	-	-	-	•	-	-	-	-	-	-	-	-	-	-	•	-	-	-	-	-	-	•	-	-	•	•	-	-	-	-	
	Programmier-Adapter 13	-	-	-	-	•	-	•	•	-	-	-	-	-	-	-	-	•	-	•	•	•	-	-	-	-	-	•	-	-	-	
	Flexstrip-Adapter	-	-	-	-	-	-	-	-	-	•	-	•	•	-	-	-	-	-	-	-	-	-	-	•	•	-	-	-	•	•	
	Programmier-Kabel CS44	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	•	-	•	-	-	-	-	-	-	-	-	-	-	-	
	Programmier-Schuh + CS44	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	•	-	-	
DRAHTLOSE PROGRAMMIERER	Noahlink Wireless	•	-	-	•	•	•	•	•	•	-	-	-	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
	ConnexxAir	-	•	-	-	-	-	-	-	-	•	•	•*	-	-	•	•	•	•	•	•	•	•	•	•	•*	-	-	-	-	-	-
	ConnexxLink	-	-	-	-	-	-	-	-	-	-	-	-	-	-	•	•	•	•	•	•	•	•	•	•*	-	-	-	-	-	-	-

\*Nur mit e2e wireless Technologie.